

Nkx2-2 Antibody, Biotin conjugated



PACO28577

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as a transcriptional activator. Required for the maintenance of NEUROD1 expression in the hormone-producing endocrine cells of the pancreas. May be involved in specifying diencephalic neuromeric boundaries, and in controlling the expression of genes that play a role in axonal guidance. Associates with chromatin at the NEUROD1 promoter region. Binds to a subset of consensus elements within the NEUROD1 promoter.

Gene ID:

Nkx2-2

Uniprot

P42586

Synonyms:

Homeobox protein Nkx-2.2 (Homeobox protein NK-2 homolog B), Nkx2-2, Nkx-2.2 Nkx2b

Immunogen:

Recombinant Mouse Homeobox protein Nkx-2.2 protein (1-273AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A