

NME6 Antibody, HRP conjugated



PACO28583

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Major role in the synthesis of nucleoside triphosphates other than ATP. The ATP gamma phosphate is transferred to the NDP beta phosphate via a ping-pong mechanism, using a phosphorylated active-site intermediate. Inhibitor of p53-induced apoptosis.

Gene ID:

NME6

Uniprot

O75414

Synonyms:

Nucleoside diphosphate kinase 6 (NDK 6) (NDP kinase 6) (EC 2.7.4.6) (Inhibitor of p53-induced apoptosis-alpha) (IPIA-alpha) (nm23-H6), NME6

Immunogen:

Recombinant Human Nucleoside diphosphate kinase 6 protein (1-186AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A