

NOL10 Antibody, Biotin conjugated



PACO28589

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleolus, poly(A) RNA binding.

Gene ID:

NOL10

Uniprot

Q9BSC4

Synonyms:

Nucleolar protein 10, NOL10

Immunogen:

Recombinant Human Nucleolar protein 10 protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A