

# NPBWR2 Antibody, HRP conjugated



PACO28607

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Interacts specifically with a number of opioid ligands. Receptor for neuropeptides B and W, which may be involved in neuroendocrine system regulation, food intake and the organization of other signals.

**Gene ID:**

NPBWR2

**Uniprot**

P48146

**Synonyms:**

Neuropeptides B/W receptor type 2 (G-protein coupled receptor 8), NPBWR2, GPR8

**Immunogen:**

Recombinant Human Neuropeptides B/W receptor type 2 protein (1-45AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A