

# NPHS1 Antibody, Biotin conjugated



PACO28621

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Seems to play a role in the development or function of the kidney glomerular filtration barrier. Regulates glomerular vascular permeability. May anchor the podocyte slit diaphragm to the actin cytoskeleton. Plays a role in skeletal muscle formation through regulation of myoblast fusion.

**Gene ID:**

NPHS1

**Uniprot**

O60500

**Synonyms:**

Nephrin (Renal glomerulus-specific cell adhesion receptor), NPHS1, NPHN

**Immunogen:**

Recombinant Human Nephrin protein (23-257AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A