

XRCC6 Antibody, Biotin conjugated



PACO28669

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

XRCC6

Uniprot

P12956

Synonyms:

X-ray repair cross-complementing protein 6 (EC 3.6.4. -) (EC 4.2.99. -) (5'-deoxyribose-5-phosphate lyase Ku70) (5'-dRP lyase Ku70) (70 kDa subunit of Ku antigen) (ATP-dependent DNA helicase 2 subunit 1) (ATP-dependent DNA helicase II 70 kDa subunit) (CTC box-binding factor 75 kDa subunit)

Immunogen:

Recombinant Human X-ray repair cross-complementing protein 6 ptotein (2-609AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A