

# NUP210 Antibody, HRP conjugated



PACO28675

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Nucleoporin essential for nuclear pore assembly and fusion, nuclear pore spacing, as well as structural integrity.

**Gene ID:**

NUP210

**Uniprot**

Q8TEM1

**Synonyms:**

Nuclear pore membrane glycoprotein 210 (Nuclear pore protein gp210) (Nuclear envelope pore membrane protein POM 210) (POM210) (Nucleoporin Nup210) (Pore membrane protein of 210 kDa), NUP210, KIAA0906

**Immunogen:**

Recombinant Human Nuclear pore membrane glycoprotein 210 protein (28-238AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A