

OGN Antibody, Biotin conjugated



PACO28713

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Induces bone formation in conjunction with TGF-beta-1 or TGF-beta-2.

Gene ID:

OGN

Uniprot

P20774

Synonyms:

Mimecan (Osteoglycin) (Osteoinductive factor) (OIF), OGN, OIF SLRR3A

Immunogen:

Recombinant Human Mimecan protein (21-298AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A