

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Odorant receptor which may be involved in sperm chemotaxis. Bourgeonal is a strong chemoattractant for sperm in vitro and is shown to be a strong agonist for OR1D2 in vitro. May also function in olfactory reception.

Gene ID:

OR1D2

Uniprot

P34982

Synonyms:

Olfactory receptor 1D2 (Olfactory receptor 17-4) (OR17-4) (Olfactory receptor OR17-6) (Olfactory receptor-like protein HGMP07E), OR1D2, OLF1

Immunogen:

Recombinant Human Olfactory receptor 1D2 protein (159-196AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A