

# OR3A1 Antibody



PACO28742

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Odorant receptor.

**Gene ID:**

OR3A1

**Uniprot**

P47881

**Synonyms:**

Olfactory receptor 3A1 (Olfactory receptor 17-40) (OR17-40) (Olfactory receptor OR17-15), OR3A1, OLFRA03

**Immunogen:**

Recombinant Human Olfactory receptor 3A1 protein (162-199AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A