

OR3A1 Antibody, HRP conjugated



PACO28743

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Odorant receptor.

Gene ID:

OR3A1

Uniprot

P47881

Synonyms:

Olfactory receptor 3A1 (Olfactory receptor 17-40) (OR17-40) (Olfactory receptor OR17-15), OR3A1, OLFRA03

Immunogen:

Recombinant Human Olfactory receptor 3A1 protein (162-199AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A