

NDUFA1 Antibody, HRP conjugated



PACO28747

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

NDUFA1 (NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 1) is an essential component of the multisubunit NADH ubiquinone oxidoreductase (complex 1), the first enzyme complex in the mitochondrial respiratory chain. Complex I transfers electrons from NADH to the respiratory chain via ubiquinone.

Gene ID:

NDUFA1

Uniprot

O15239

Synonyms:

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 (Complex I-MWFE) (CI-MWFE) (NADH-ubiquinone oxidoreductase MWFE subunit), NDUFA1

Immunogen:

Recombinant Human NADH dehydrogenase [ubiquinone] 1 α subcomplex subunit 1 protein (1-70AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A