

# PAFAH1B2 Antibody, HRP conjugated



PACO28791

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Inactivates PAF by removing the acetyl group at the sn-2 position. This is a catalytic subunit.

**Gene ID:**

PAFAH1B2

**Uniprot**

P68402

**Synonyms:**

Platelet-activating factor acetylhydrolase IB subunit beta (EC 3.1.1.47) (PAF acetylhydrolase 30 kDa subunit) (PAF-AH 30 kDa subunit) (PAF-AH subunit beta) (PAFAH subunit beta), PAFAH1B2, PAFAHB

**Immunogen:**

Recombinant Human Platelet-activating factor acetylhydrolase IB subunit &beta; protein (2-229AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A