

PANK1 Antibody



PACO28798

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Plays a role in the physiological regulation of the intracellular CoA concentration.

Gene ID:

PANK1

Uniprot

Q8TE04

Synonyms:

Pantothenate kinase 1 (hPank) (hPank1) (EC 2.7.1.33) (Pantothenic acid, kinase 1),
PANK1, PANK


Immunogen:

Recombinant Human Pantothenate kinase 1 protein (462-582AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO28798 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO28798 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).