

PAPLN Antibody



PACO28802

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PAPLN

Uniprot

O95428

Synonyms:

Papilin, PAPLN

Immunogen:

Recombinant Human Papilin protein (19-129AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A