

PCBP4 Antibody



PACO28830

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Single-stranded nucleic acid, binding protein that binds preferentially to oligo dC.

Gene ID:

PCBP4

Uniprot

P57723

Synonyms:

Poly(rC)-binding protein 4 (Alpha-CP4), PCBP4

Immunogen:

Recombinant Human Poly(rC)-binding protein 4 protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A