

PCOLCE2 Antibody, FITC conjugated



PACO28840

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to the C-terminal propeptide of types I and II procollagens and may enhance the cleavage of that propeptide by BMP1.

Gene ID:

PCOLCE2

Uniprot

Q9UKZ9

Synonyms:

Procollagen C-endopeptidase enhancer 2 (Procollagen COOH-terminal proteinase enhancer 2) (PCPE-2) (Procollagen C-proteinase enhancer 2), PCOLCE2, PCPE2

Immunogen:

Recombinant Human Procollagen C-endopeptidase enhancer 2 protein (24-415AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A