

# PCP4 Antibody



PACO28842

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable regulator of calmodulin signaling.

**Gene ID:**

PCP4

**Uniprot**

P48539

**Synonyms:**

Calmodulin regulator protein PCP4 (Brain-specific polypeptide PEP-19) (Purkinje cell protein 4), PCP4, PEP19

**Immunogen:**

Recombinant Human Calmodulin regulator protein PCP4 protein (1-62AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A