

# PFDN2 Antibody, Biotin conjugated



PACO28893

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds specifically to cytosolic chaperonin (c-CPN) and transfers target proteins to it. Binds to nascent polypeptide chain and promotes folding in an environment in which there are many competing pathways for nonnative proteins.

**Gene ID:**

PFDN2

**Uniprot**

Q9UHV9

**Synonyms:**

Prefoldin subunit 2, PFDN2, PFD2

**Immunogen:**

Recombinant Human Prefoldin subunit 2 protein (1-154AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A