

# Plaur Antibody, Biotin conjugated



PACO28989

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA.

**Gene ID:**

Plaur

**Uniprot**

P49616

**Synonyms:**

Urokinase plasminogen activator surface receptor (U-PAR) (uPAR) (CD\_antigen: CD87), Plaur

**Immunogen:**

Recombinant Rat Urokinase plasminogen activator surface receptor ptotein (25-299AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A