

# PNMA1 Antibody, Biotin conjugated



PACO29037

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, focal adhesion, nucleolus, inflammatory response to antigenic stimulus, positive regulation of apoptotic process.

**Gene ID:**

PNMA1

**Uniprot**

Q8ND90

**Synonyms:**

Paraneoplastic antigen Ma1 (37 kDa neuronal protein) (Neuron- and testis-specific protein 1), PNMA1, MA1

**Immunogen:**

Recombinant Human Paraneoplastic antigen Ma1 protein (151-353AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A