

POP4 Antibody, HRP conjugated



PACO29093

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Part of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends. May function with RPP38 to coordinate the nucleolar targeting and/or assembly of RNase P.

Gene ID:

POP4

Uniprot

O95707

Synonyms:

Ribonuclease P protein subunit p29 (hPOP4) (EC 3.1.26.5), POP4, RPP29

Immunogen:

Recombinant Human Ribonuclease P protein subunit p29 protein (1-220AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A