

PRAMEF5 Antibody, Biotin conjugated



PACO29155

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	PRAMEF5
Source:	Uniprot
Rabbit	Q5TYX0
Isotype:	Synonyms:
IgG	PRAME family member 5, PRAMEF5, PRAMEF23 PRAMEF5L
Applications:	Immunogen:
ELISA	Recombinant Human PRAME family member 5 protein (1-476AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A