

PACO29196

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Cytoplasm, extracellular exosome, intracellular, pyridoxal phosphate binding.

**Gene ID:**

PLPBP

**Uniprot**

O94903

**Synonyms:**

Pyridoxal phosphate homeostasis protein (PLP homeostasis protein) (Proline synthase co-transcribed bacterial homolog protein) (Pyridoxal phosphate-binding protein), PLPBP, PROSC

**Immunogen:**


Recombinant Human Pyridoxal phosphate homeostasis protein (1-275AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO29196 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO29196 at dilution of 1:100.