

RAB27B Antibody



PACO29320

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in targeting uroplakins to urothelial apical membranes.

Gene ID:

RAB27B

Uniprot

O00194

Synonyms:

Ras-related protein Rab-27B (C25KG), RAB27B

Immunogen:

Recombinant Human Ras-related protein Rab-27B protein (2-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A