

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in melanosomal transport and docking. Involved in the proper sorting of TYRP1. Involved in peripheral melanosomal distribution of TYRP1 in melanocytes; the function, which probably is implicating vesicle-trafficking, includes cooperation with ANKRD27 and VAMP7. Plays a role in the maturation of phagosomes that engulf pathogens, such as S.aureus and M. tuberculosis.

Gene ID:

RAB38

Uniprot

P57729

Synonyms:

Ras-related protein Rab-38 (Melanoma antigen NY-MEL-1), RAB38

Immunogen:

Recombinant Human Ras-related protein Rab-38 protein (1-211AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A