

RCOR3 Antibody



PACO29428

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May act as a component of a corepressor complex that represses transcription.

Gene ID:

RCOR3

Uniprot

Q9P2K3

Synonyms:

REST corepressor 3, RCOR3, KIAA1343


Immunogen:

Recombinant Human REST corepressor 3 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO29428 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO29428 at dilution of 1:100.