

RGS22 Antibody



PACO29488

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits thereby driving them into their inactive GDP-bound form.

Gene ID:

RGS22

Uniprot

Q8NE09

Synonyms:

Regulator of G-protein signaling 22 (RGS22), RGS22

Immunogen:

Recombinant Human Regulator of G-protein signaling 22 protein (201-504AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A