

S100g Antibody



PACO29728

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

S100g

Uniprot

P02634

Synonyms:

Protein S100-G (9 kDa CaBP) (Calbindin-D9k) (Cholecalciferol) (S100 calcium-binding protein G) (Vitamin D-dependent calcium-binding protein, intestinal) (CABP), S100g, Calb3 S100d

Immunogen:

Recombinant Rat Protein S100-G protein (2-79AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A