

SCNM1 Antibody



PACO29764

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a role in RNA splicing, possibly contributing to the recognition of non-consensus donor sites.

Gene ID:

SCNM1

Uniprot

Q9BWG6

Synonyms:

Sodium channel modifier 1, SCNM1


Immunogen:

Recombinant Human Sodium channel modifier 1 protein (1-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO29764 at dilution of 1:100.