

# Serpina1 Antibody, FITC conjugated



PACO29830

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Inhibitor of serine proteases. The primary target is elastase, but also has a moderate affinity for plasmin and thrombin.

**Gene ID:**

Serpina1

**Uniprot**

P17475

**Synonyms:**

Alpha-1-antiproteinase (Alpha-1-antitrypsin) (Alpha-1-proteinase inhibitor) (Serpina A1), Serpina1

**Immunogen:**

Recombinant Rat Alpha-1-antiproteinase protein (25-411AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A