

Serpinf1 Antibody, Biotin conjugated



PACO29851

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Serpinf1

Uniprot

Q80ZA3

Synonyms:

Pigment epithelium-derived factor (PEDF) (Caspin) (Serpin F1) (Stromal cell-derived factor 3) (SDF-3), Serpinf1, Pedf Sdf3

Immunogen:

Recombinant Rat Alpha-2 antiplasmin protein (20-418AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A