

# HLA-G Antibody, HRP conjugated



PACO29857

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the presentation of foreign antigens to the immune system. Plays a role in maternal tolerance of the fetus by mediating protection from the deleterious effects of natural killer cells, cytotoxic T-lymphocytes, macrophages and mononuclear cells.

**Gene ID:**

HLA-G

**Uniprot**

P17693

**Synonyms:**

HLA class I histocompatibility antigen, alpha chain G (HLA G antigen) (MHC class I antigen G), HLA-G, HLA-6.0 HLAG

**Immunogen:**

Recombinant Human HLA class I histocompatibility antigen, &alpha; chain G protein (25-337AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A