

SFTA3 Antibody



PACO29864

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SFTA3

Uniprot

P0C7M3

Synonyms:

Surfactant-associated protein 3 (Surfactant-associated protein H) (SP-H), SFTA3, SFTPH

Immunogen:

Recombinant Human Surfactant-associated protein 3 protein (1-94AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A