

# SH3BGR Antibody



PACO29876

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytosol, SH3/SH2 adaptor activity, positive regulation of signal transduction, protein complex assembly.

**Gene ID:**

SH3BGR

**Uniprot**

P55822

**Synonyms:**

SH3 domain-binding glutamic acid, rich protein (SH3BGR protein) (21-glutamic acid, rich protein) (21-GARP), SH3BGR

**Immunogen:**

Recombinant Human SH3 domain-binding glutamic acid, rich protein (1-239AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A