

SH3BGRL Antibody



PACO29880

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, extracellular exosome, extracellular space, nucleus, SH3/SH2 adaptor activity, positive regulation of signal transduction.

Gene ID:

SH3BGRL

Uniprot

O75368

Synonyms:

SH3 domain-binding glutamic acid, rich-like protein, SH3BGRL


Immunogen:

Recombinant Human SH3 domain-binding glutamic acid, rich-like protein (1-114AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO29880 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO29880 at dilution of 1:100.