

SLC31A2 Antibody, Biotin conjugated



PACO29951

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in low-affinity copper uptake.

Gene ID:

SLC31A2

Uniprot

O15432

Synonyms:

Probable low affinity copper uptake protein 2 (Copper transporter 2) (hCTR2) (Solute carrier family 31 member 2), SLC31A2, COPT2 CTR2

Immunogen:

Recombinant Human Probable low affinity copper uptake protein 2 protein (44-93AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A