

# SMCP Antibody, Biotin conjugated



PACO30003

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in sperm motility. Its absence is associated with genetic background dependent male infertility. Infertility may be due to reduced sperm motility in the female reproductive tract and inability to penetrate the oocyte zona pellucida.

**Gene ID:**

SMCP

**Uniprot**

P49901

**Synonyms:**

Sperm mitochondrial-associated cysteine-rich protein, SMCP, MCS MCSP

**Immunogen:**

Recombinant Human Sperm mitochondrial-associated cysteine-rich protein (1-116AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A