

SNCG Antibody, Biotin conjugated



PACO30027

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in neurofilament network integrity. May be involved in modulating axonal architecture during development and in the adult. In vitro, increases the susceptibility of neurofilament-H to calcium-dependent proteases By similarity. May also function in modulating the keratin network in skin. Activates the MAPK and Elk-1 signal transduction pathway By similarity.

Gene ID:

SNCG

Uniprot

O76070

Synonyms:

Gamma-synuclein (Breast cancer-specific gene 1 protein) (Persyn) (Synoretin) (SR), SNCG, BCSG1 PERSYN PRSN

Immunogen:

Recombinant Human Gamma-synuclein protein (1-127AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A