

SNX24 Antibody



PACO30052

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in several stages of intracellular trafficking.

Gene ID:

SNX24

Uniprot

Q9Y343

Synonyms:

Sorting nexin-24, SNX24


Immunogen:

Recombinant Human Sorting nexin-24 protein (1-169AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO30052 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO30052 at dilution of 1:100.