

# SOD3 Antibody, HRP conjugated



PACO30081

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Protect the extracellular space from toxic effect of reactive oxygen intermediates by converting superoxide radicals into hydrogen peroxide and oxygen.

**Gene ID:**

SOD3

**Uniprot**

P08294

**Synonyms:**

Extracellular superoxide dismutase [Cu-Zn] (EC-SOD) (EC 1.15.1.1), SOD3

**Immunogen:**

Recombinant Human Extracellular superoxide dismutase [Cu-Zn] protein (19-160AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A