

# SPRY1 Antibody, FITC conjugated



PACO30118

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May function as an antagonist of fibroblast growth factor (FGF) pathways and may negatively modulate respiratory organogenesis.

**Gene ID:**

SPRY1

**Uniprot**

O43609

**Synonyms:**

Protein sprouty homolog 1 (Spry-1), SPRY1

**Immunogen:**

Recombinant Human Protein sprouty homolog 1 protein (1-178AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A