

HPCAL1 Antibody, Biotin conjugated



PACO30143

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in the calcium-dependent regulation of rhodopsin phosphorylation.

Gene ID:

HPCAL1

Uniprot

P37235

Synonyms:

Hippocalcin-like protein 1 (Calcium-binding protein BDR-1) (HLP2) (Visinin-like protein 3) (VILIP-3), HPCAL1, BDR1

Immunogen:

Recombinant Human Hippocalcin-like protein 1 protein (2-193AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A