

SSB Antibody, Biotin conjugated



PACO30147

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds to the 3' poly(U) terminus of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and facilitating their folding and maturation.

Gene ID:

SSB

Uniprot

P05455

Synonyms:

Lupus La protein (La autoantigen) (La ribonucleoprotein) (Sjogren syndrome type B antigen) (SS-B), SSB

Immunogen:

Recombinant Human Lupus La protein (1-408AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A