

STX10 Antibody, FITC conjugated



PACO30200

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network.

Gene ID:

STX10

Uniprot

O60499

Synonyms:

Syntaxin-10 (Syn10), STX10, SYN10

Immunogen:

Recombinant Human Syntaxin-10 protein (1-179AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A