

STX6 Antibody, Biotin conjugated



PACO30205

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in intracellular vesicle trafficking.

Gene ID:

STX6

Uniprot

O43752

Synonyms:

Syntaxin-6, STX6

Immunogen:

Recombinant Human Syntaxin-6 protein (5-234AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A