

SUCNR1 Antibody, HRP conjugated



PACO30215

Product Information

Size:	Protein Background:
50ug	Receptor for succinate.
Reactivity:	Gene ID:
Human	SUCNR1
Source:	Uniprot
Rabbit	Q9BXA5
Isotype:	Synonyms:
IgG	Succinate receptor 1 (G-protein coupled receptor 91) (P2Y purinoceptor 1-like), SUCNR1, GPR91
Applications:	Immunogen:
ELISA	Recombinant Human Succinate receptor 1 protein (303-334AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A