

SUPT3H Antibody, FITC conjugated



PACO30232

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcriptional activator.

Gene ID:

SUPT3H

Uniprot

O75486

Synonyms:

Transcription initiation protein SPT3 homolog (SPT3-like protein), SUPT3H, SPT3

Immunogen:

Recombinant Human Transcription initiation protein SPT3 homolog protein (1-328AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A