

# SYF2 Antibody, HRP conjugated



PACO30247

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	May be involved in pre-mRNA splicing.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	SYF2
<b>Source:</b>	<b>Uniprot</b>
Rabbit	O95926
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Pre-mRNA-splicing factor SYF2 (CCNDBP1-interactor) (p29), SYF2, CBPIN GCIPIP
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Pre-mRNA-splicing factor SYF2 protein (2-243AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A