

# TBCA Antibody, Biotin conjugated



PACO30285

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Tubulin-folding protein; involved in the early step of the tubulin folding pathway.

**Gene ID:**

TBCA

**Uniprot**

O75347

**Synonyms:**

Tubulin-specific chaperone A (TCP1-chaperonin cofactor A) (Tubulin-folding cofactor A) (CFA), TBCA

**Immunogen:**

Recombinant Human Tubulin-specific chaperone A protein (2-108AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A